


1

## Objectives


1. What is a visual field?
2. What is visual field testing?
3. How to perform a visual field test
4. What are the different types of visual fields?
5. What are different types of visual field defects?
6. How to get good results!



2



## What is a visual field?

Visual field is the area that can be seen without moving the eyes or head<sup>23,45</sup>. It includes what can be seen straight ahead and to the sides<sup>24</sup>. Vision is usually best in the middle of the visual field<sup>24</sup>. A visual field test can measure the extent and quality of the visual field as part of an eye exam<sup>14</sup>.




3

## Your Visual Field




4



5

## How many cranial nerves are there

- How many total cranial nerves are there?




6


The Neuron System is complex

- 12 Cranial nerves
- 6 of them have a significant impact on the ocular system
- CN #2 is the optic nerve
- Vision occurs in the occipital lobe of the brain


How parts of the brain got their names: ...see more



7




8



Visual Pathway

- Visual pathway has seven structures
- Retina
- Optic Nerve – light super-highway information
- Optic Chiasm
- Optic Tract
- Lateral Geniculate Body (LGB)
- Optic Radiations
- Visual Cortex ...where vision occurs

9



What is visual field testing?

- As you focus on the words in this article, how much can you see out of the corners of your eyes? Can you tell what's happening in your surroundings?
- **Your visual field is how wide of an area your eye can see when you focus on a central point.**
- Visual field testing is one way your Eye Care Provider measures how much vision you have in either eye, and how much vision loss may have occurred over time.

10

Two types of Visual Field Tests

- Kinetic
  - Uses moving light targets
  - Goldmann Field Perimetry
- Static
  - Automated with blinking lights
  - Brightness (0 = Brightest and 50 = dimmest)
  - Different target sizes (1-5)
  - Different speeds (normal or slow)
  - Different programs
    - 10-2 for central field 68 points (Plaqueinil & vein/artery occlusions)
    - 24-2 standard glaucoma 54 points
    - 30-2 neuro pathway testing 76 points
    - 36 point – Ptosis testing
    - 120 point – screener



11

Glaucoma Testing Terms

- ☐ DFE
- ☐ FDT
- ☐ Heidelberg (new technology)
- ☐ OCT
- ☐ HVF (reduces operator error)
- ☐ Threshold testing
- ☐ FastPac Testing
- ☐ Rx is important
- ☐ Pupil Diameter (3mm min)
- ☐ Numeric data
- ☐ Gray scale
- ☐ Total deviation
- ☐ Pattern deviation
- ☐ Global Hemifield Test (GHT)
- ☐ Global indices (standards)
- ☐ Pattern Standard Deviation (PSD)
- ☐ Mean Deviation (MD)

12


### Different types of visual field testing?

- Confrontation fields
- Automatic
  - 10-2 Macula, 24-2 Glaucoma, 30-2 – Neuro, 120 pt Screener, 36 pt Ptosis,
- Amsler Grid (key chain version)
- Static
- Kinetic
- FDT – Frequency Doubling Technology
- Headset Mounted




13

### Automated Visual Field – Industry Standard



14



How to perform Visual Field Testing

- Complete pt demographic info
- Add SRx and calculate age/Rx for test
- Select the appropriate test
- Patch the left eye\*
- Proper pt instructions
- Proper alignment (pupil size 3mm)
- Monocular Field
- Add SRx and calculate \*age
- Begin
- Different fields (SNIT)
  - Superior - 60 degrees
  - Nasal - 60 degrees
  - Inferior - 70 degrees
  - Temporal - 100 degrees

15

### Visual Field Testing

- Enter patient data
  - Recall patient data for establish pt
  - Enter patient's spectacle Rx
  - Add the correction for OD in lens holder\*
- Select the appropriate test
- Position patient, check for comfort (**critical for test results**)
  - Proper alignment with monitor crosshair
- Constant monitoring (**never leave the patient alone**)
- Provide the appropriate instructions
- Complete test – save results!



16

### Visual Field Testing

- If the patient is dilated add +3.00 for lack of accommodation
- If the patient has had cataract surgery with standard lenses add +3.25
- When unable to perform test document the chart to reflect



17



### Patient Alignment and Instructions

- This test will determine the sensitivity of the connection between the eye and the brain
- The test will measure your central and peripheral vision
- You are to maintain your fixation on the fixation target
- There will be lights around the fixation target, do not chase them with your eyes

18

Patient instructions continued...


- Some of the lights will be bright, some dim barely visible, when ever you see a light flash, hit the responder
- If you need a break, please hold the button down and do not let it go until we talk
- Blink normally while testing
- If you get more than 3 fixation losses, we will restart the test



19

Analysis of Automated Visual Field Tests


- Reliability is important
- Threshold Values
- Grayscale map – patient sensitivity
- Total deviation – depth of field
- Pattern deviation – localized defects
- Glaucoma Hemifield Test (GHT) above and below the horizontal line
- Visual Field Indices -
- Probability Symbols
- Gaze Tracking



20

Analysis of Automated Visual Field Tests


- SITA Fast or SITA Standard
- Reliability is important
- Threshold Values
- Grayscale map – pictorial overview of PSE raw data
- Total deviation decibels based on age
- Pattern deviation – factors out relative noise, sine waves
- Glaucoma Hemifield Test (GHT) – takes advantage of asymmetry field loss
- Visual Field Indices
- Probability Symbols
- Gaze Tracking



21

Interpretation Guide


- Reliability Indices (FP, FN, FL)
- **False positives** (trigger happy)
    - Responded when there is no stimulus
  - **False negatives**
    - Missing previous area responded to
  - **Fixation loss rate\***
    - Patient not fixated on the target
    - Poor blinking\*
  - Gray scale map
    - Raw decibel sensitivity / dark area –reduced sensitivity
  - Threshold Values
    - Measured decibel – test point sensitivity




22

Interpretation Guide cont...

- Total and Pattern deviation
  - Must compare total and pattern deviation maps to each other
    - Numerical total deviation
    - Numerical pattern deviation
- Total probability -
- Visual Field Indices
  - VFI – less affected by cataracts / ganglion cell layer
  - MD - age-normal deviation
  - PSD – irregularities in the field from localized field loss
- Glaucoma Hemifield Test (GHT)
  - Pattern deviation associated with glaucoma



23



24

### Goldmann Perimetry

In Goldmann perimetry, an isopter is plotted kinetically by moving the test target from the periphery or non-seeing area to an area where it is seen by the patient. The test target is also projected statically at a single location and the brightness increased until the patient sees it. These two methods are then used to determine a patient's field of vision, or the actual area within which the eye is able to detect the presence of an object or stimulus.



25

**Reliability Indices (FP, FN, FL)**


- False positives** (trigger happy)
  - Responded when there is no stimulus
- False negatives**
  - Missing previous response area seen
  - PT may be dozing off
- Fixation loss rate**
  - Patient not fixated on the target
  - In the blind spot area
- Gray scale map**
  - Raw decibel sensitivity / dark area -reduced sensitivity
- Threshold Values**
  - Measured decibel - test point sensitivity
  - 0-50 0=brightest 50=dimmiest

### Interpretation Guide

26

- Total and Pattern deviation**
  - Must compare total and pattern deviation maps to each other
    - Numerical total deviation
    - Numerical pattern deviation
  - Compares to normative data base
- Total probability -**
- Visual Field Indices**
  - VFI - less affected by cataracts / ganglion cell layer
  - MD - age-normal deviation
  - PSD - irregularities in the field from localized field loss
- Glaucoma Hemifield Test (GHT)**
  - Pattern deviation associated with glaucoma

### Interpretation Guide cont...



27

### Virtual Visual Field


- Consistent with virtual reality instructions
- Better with patients with tremors
- Better for patients with back and neck issues
- Caution with claustrophobic pts
- Caution with patients with vertigo
- We can't see the eyes of the pts



28

### Goldmann Perimetry

In Goldmann perimetry, an isopter is plotted kinetically by moving the test target from the periphery or non-seeing area to an area where it is seen by the patient. The test target is also projected statically at a single location and the brightness increased until the patient sees it. These two methods are then used to determine a patient's field of vision, or the actual area within which the eye can detect the presence of an object or stimulus.



29


### 3D Headsets

Meta Quest 2




30

### Different types of visual field defects?



31


### Different types of visual field defects?




32

### Visual Field Defects

- Common types of field defects
- Blind spots - Areas of blindness in the visual field
- Hemianopsia - Blindness in one half of the visual field of one or both eyes
  - **Homonymous** - Involving the **nasal** half of the visual field of one eye and the **temporal** half of the visual field of the other eye
  - **Incongruous** - strange shaped visual fields




33



34

### Types of Visual Fields


- 10-2 - central field test
- 24-2 - std glaucoma test
- 30-2 - std neurological test
- 36pt - ptosis test
- 120 pt - screening field
- Drivers license exam



35


### Conditions where VF testing useful

- MS
- Stroke
- Optic Neuritis
- Retinitis Pigmentosa
- Vein/Arterial occlusion
- Plaque toxicity



36

### Emphasis On The Macula – OCT Testing



- Central 10-2 for early-stage glaucoma
- Central 12 points of the 24-2 does cover most patients

37

### When test is complete

- Have the patient sit back
- Verify testing is within the protocol parameters
- Ensure test is saved
- Return the lenses to the trial lens set
- **Verify entry into pts medical chart**
- Bill the testing visit



38

### Making Sure The You Get Good VF Results

- Staff training is a great place to begin, don't assume they know
  - Was the test done on SITA Standard or SITA Fast? (Swedish Interactive Threshold Algorithm)
- Verify your patient is the correct patient
- Pupil diameter must be at least 3mm
- Verify the correct test (don't assume the provider is ordering the same tests)
- Patient position is paramount, comfort is king!
- Always monitor, never leave them alone, encourage them thru testing
  - Tell them if there is a fixation loss, and tell them if you have to restart test
- You may need to change test speed if patient response time is slow
- If patient has DED you may need to give lubricating drops before start
  - If the patient is tearing, you may need to pause the testing
- You may need to pause the testing if you notice pt struggling
- If the patient is wiggling, find out if they are comfortable

39

### When not to perform a visual field test!


- When the patient is not well!
- If the patient says they are in too much pain to test
- Visual acuity is:
  - Hand motion (HM)
  - Light perception (LP)
  - No light perception (NLP)
- If there are no orders for the test in the patient's chart

40

### Billing examples

- Plaquenil
- Glaucoma
- Neuro testing
- Ptoisis testing

41



[martralyn@msn.com](mailto:martralyn@msn.com)

42